


# Ancillary Reagents

Simplify & Standardize Your FISH Workflows with Empire Ancillary Reagents



## Complete Line of Ancillary Reagents

Empire Genomics is excited to introduce our latest ancillary reagents, designed to enhance your genomic research. Our new lineup includes DAPI with Antifade for superior fluorescence preservation, Hot & Cold Wash Buffers for optimized washing protocols, FFPE Pre-Treatment and Cyto Pre-Treatment Kits for efficient sample preparation and both SwiftFISH and standard Hybridization Buffers for streamlined hybridization processes. These innovative products are crafted to provide reliability and exceptional performance in your laboratory workflows.

- ▶ DAPI w/ Antifade
- ▶ Hot & Cold Wash Buffers
- ▶ FFPE Pre-Treatment Kit
- ▶ Cyto Pre-Treatment Kit
- ▶ Rapid & Standard Hybridization Buffers

**Distributed by:**

**CliniSciences Group**



## Pre-Treatment Kits

Proper pretreatment is crucial for the success of Fluorescence In Situ Hybridization (FISH) assays. Our FFPE and Cyto Pretreatment Kits ensure optimal sample preparation by enhancing probe accessibility and improving hybridization efficiency. These steps are vital for achieving accurate, reliable results in genomic mapping, cancer diagnostics, and genetic disease detection. Trust our pretreatment solutions to deliver the clarity and precision your FISH assays demand.

### FFPE Pre-Treatment Kit

The FFPE Pre-Treatment Kit is designed for optimal heat pretreatment and digestion of formalin-fixed, paraffin-embedded (FFPE) tissue slides. This comprehensive kit simplifies sample preparation, ensuring effective tissue digestion and high-quality results for your genomic research and Laboratory Developed Tests (LDT).

Product Details		
Application & Volume	Tissue Based FISH and CISH, 25 Tests	
Kit Components	Size	Concentration
Pretreatment Solution	250 ml	10mM Citric Acid
Wash Buffer	250 ml	2 x SSC
Protease Buffer	250 ml	0.01N HCl
Protease	5 x 20mg	Pepsin solid

### Cyto Pre-Treatment Kit

Empire Genomics Cyto Pre-Treatment Kit is designed for the efficient digestion of thin prep urine slides, ensuring optimal sample quality for cytogenetic analysis. This kit simplifies the preparation process, delivering consistent and high-quality results for your genomic research and Laboratory Developed Tests (LDT).

Product Details		
Application & Volume	Urine Based FISH and CISH, 25 Tests	
Kit Components	Size	Concentration
Wash Buffer	250 ml	2 x SSC
Protease Buffer	250 ml	0.01N HCl
Protease	5 x 40mg	Pepsin solid

**Distributed by:**

**CliniSciences Group**

## Buffers

Buffers play a critical role in the success of Fluorescence In Situ Hybridization (FISH) assays by ensuring proper probe hybridization and signal clarity. Our comprehensive range includes Hot and Cold Wash Buffers for thorough sample washing, Rapid Hybridization Buffer (SwiftFISH) for expedited hybridization processes, and Standard Hybridization Buffer for consistent and reliable results. These buffers are designed to enhance the performance of your FISH assays, providing the precision and accuracy needed for advanced genomic research and diagnostics.

### Hot & Cold Wash Buffers **Top Seller**

Our Wash Buffer I (0.3% Igepal/0.4x SSC) and Wash Buffer II (0.1% Igepal/2x SSC) are sterile solutions intended for use as a post-hybridization hot and cold wash. Suitable for washing both FISH and CISH slides, this product is shipped in a cubitainer and includes a spout for easy dispensing.

Product Details	
Volume	4L each

### SwiftFISH Hybridization Buffer **Two Hour Turn-Around, Excellent Fluidity**

Reduce your hybridization time from an overnight cycle to as little as two hours or run it overnight while maintaining the same level of precision and signal strength, with Empire Genomics SwiftFISH® rapid hybridization buffer. The SwiftFISH buffer has been tested on FFPE, Peripheral Blood, and Bone Marrow samples and is optimized to reduce signal background.

Product Details	
Volume	200ul / 20 Tests
Application	FISH and CISH

### Hybridization Buffer

Our Hybridization Buffer is engineered to ensure optimal binding conditions for your FISH and CISH assays, delivering consistent and reliable results. This high-quality buffer enhances probe hybridization efficiency, making it an essential component.

Product Details	
Volume	200ul / 20 Tests
Application	FISH and CISH

### DAPI w/ Antifade Kit **Tech Preferred**

Experience the next level of convenience and reliability in your fluorescence in situ hybridization (FISH) experiments with our Ready-to-Use DAPI with Antifade Solution. Designed to save you time and enhance your results, this solution is ready straight from the bottle, eliminating the need for any preparation steps.

Product Details	
Volume	500ul / 50 Tests
Concentration	125ng/ml



# Ordering Information

Product	Category	SKU	Intended Use	Market
DAPI w/ Antifade	Counterstain	DAPI-0125	RUO	Global
Wash Buffer 1	Post Hybridization Hot Wash	WASH1-004L	RUO	Global
Wash Buffer 2	Post Hybridization Cold Wash	WASH2-004L	RUO	Global
FFPE Pre-Treatment Kit	Pre-Treatment	FPTK-001	RUO	Global
Cyto Pre-Treatment Kit	Pre-Treatment	CPTK-002	RUO	Global
SwiftFISH Buffer	Rapid Hybridization Buffer	SwiftFISH-Buffer	RUO	Global
Hybridization Buffer	Hybridization Buffer	Hyb-Buffer	RUO / CE	Global



For In Vitro Use Only | For Research Use Only | Not For Diagnostic Use

**Distributed by:**

**CliniSciences Group**

# CliniSciences Group

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oesterreich@clinisciences.com](mailto:oesterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [denmark@clinisciences.com](mailto:denmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Biotrend Chemikalien GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [info@biotrend.com](mailto:info@biotrend.com)  
Web: <https://www.biotrend.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Kraaijenhoffstraat 137A,  
1018RG Amsterdam, Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.Z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 91 269 40 65  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences Limited  
Address: Marktgasse 18 8302 Kloten -  
Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: Biotrend Chemicals LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [info@biotrend-usa.com](mailto:info@biotrend-usa.com)  
Web: <https://www.biotrend-usa.com>

