

KAC^TUS

MaxNuclease, GMP-Grade

DMF #036799

All-Purpose Nuclease



Distributed by:

CliniSciences Group

MaxNuclease™ for Nucleic Acid Degradation

MaxNuclease endonuclease identified from *Serratia marcescens* is genetically engineered and expressed in *E. coli* under cGMP manufacturing standards. MaxNuclease is a non-specific nuclease with high activity and specificity that degrades all forms of nucleic acids including single- and double-stranded, linear and circular nucleic acids. The enzyme is a homodimer of two 30 kDa subunits containing two disulfide bonds that are essential for activity and stability. It hydrolyzes internal phosphodiester

bonds between nucleotides in nucleic acids to produce 5'-monophosphate oligonucleotides of 2-5 bases in length. MaxNuclease is ideal for purification of viral vectors and viral vaccines, as well as for protein purification and other applications where removal of contaminating nucleic acids is desired. It can also effectively reduce the viscosity of cell lysates and limit cell aggregation and clumping.

Thorough & Efficient Removal of Nucleic Acids




Figure 1. The virus harvest solution was treated with 25U/mL and 50U/mL endonuclease at 37°C for 2h, respectively. Detection of HCD residue (left) and pDNA residue (right) was analyzed. KACTUS MaxNuclease has higher degradation activity versus Competitor B in both HCD residue and pDNA residue testing for both 25U/mL and 50U/mL working concentrations.

Ordering Information

Catalog Number	Product Name	Quantities
GMP-NUC-SE101	MaxNuclease, GMP-Grade	250KU / 5MU
NUC-SE00B	MaxNuclease ELISA Kit	96 Tests

Distributed by:

CliniSciences Group


Features

- Manufactured in a GMP-compliant facility
- Raw materials free from animal-derived components
- Strict quality management to meet clinical manufacturing standards
- FDA Drug Master Files Type II filing (DMF #036799)
- A complete document package to support your project registration

Applications

- Purification of viral vaccines and viral vectors (lentivirus, adenovirus, oncolytic virus, etc.)
- Removal of nucleic acid residues (DNA/RNA) in biological products
- Reducing the viscosity of cell lysates and cell supernatants
- Preparing samples in western blot, 2D gel electrophoresis, ELISA, and chromatography to improve resolution and recovery

Degradation of PCR Product, Genomic DNA, and Plasmid DNA



- Lane M DNA marker
- Lane 1 PCR product
- Lane 2 PCR product +1U MaxNuclease
- Lane 3 PCR product +1U competitor
- Lane 4 genomic DNA
- Lane 5 genomic DNA +1U MaxNuclease
- Lane 6 genomic DNA +1U competitor
- Lane 7 plasmid DNA
- Lane 8 plasmid DNA +1U MaxNuclease
- Lane 9 plasmid DNA +1U competitor


Figure 1. MaxNuclease can degrade any form of nucleic acid such as PCR products, gDNA, plasmids, and RNA.

MaxNuclease ELISA Kit

MaxNuclease ELISA Kit can detect and quantitatively analyze the MaxNuclease residues in viral vectors and viral vaccines with high sensitivity and specificity.

The kit uses sandwich ELISA to determine the concentration of MaxNuclease in the test sample.

Example Standard Curve



Detection Range	Sensitivity	Precision
46.88 pg/mL - 3000.00 pg/mL	23.44 pg/mL	CV < 10%

Distributed by:

CliniSciences Group

Customizable GMP Documentation Package

- Datasheet
- CoA
- CoO
- MSDS
- Melamine Statement
- TSE/BSE Statement
- Nitrosamine Statement
- DMF Filing

Product Specifications

Parameter	Specification
Source	<i>E. coli</i> with endonuclease gene from <i>Serratia marcescens</i>
Molecular Weight	Approximately 27.8 kDa
Purity	>99% by SEC-HPLC
Activity	≥250 U/μL, tested by degradation of Salmon Sperm DNA
Formulation	20mM Tris-HCl, 20mM NaCl, 2mM MgCl ₂ , 50% Glycerol, pH 8.0
Endotoxin	Less than 0.01EU/kU, determined by LAL method
Sterility	Negative
Mycoplasma	Negative, tested by qPCR
Storage	Store at -20±5°C. Avoid repeated freeze-thaw.
Unit Definition	One unit corresponds to the amount of enzyme required to produce a change in absorbance at 260 nm of 1.0 in 30 minutes, at 37°C and pH 8.0.

Reaction Conditions

Condition	Optimal*	Effective**
Mg ²⁺	1-2mM	1-10mM
Na ⁺ , K ⁺	0-100mM	0-300mM
pH	8-10	4-10
Temperature	37°C	0-50°C
PO ₄ ³⁻	0-10mM	0-80mM

*Optimal is defined as the condition under which MaxNuclease retains >90% of its activity

**Effective is defined as the condition under which MaxNuclease retains >15% of its activity

Distributed by:

CliniSciences Group

Effect of cations on enzyme activity




Figure 2. A minimum concentration of 1 to 2 mM Mg²⁺ is essential for activity of MaxNuclease.




Figure 3. Na⁺ and K⁺ strongly inhibit the endonuclease activity. Activity is lost when the concentrations reach 500 mM.

Effect of Reaction Temperature and pH on Enzyme Activity




Figure 4. Effect of temperature on MaxNuclease endonuclease activity. The relative activity rises with increasing temperature. The optimum temperature is 37 °C.




Figure 5. Effect of pH on MaxNuclease endonuclease activity. The optimum pH is between 8 and 10.

Effect of Common Buffers on Enzyme Activity




Figure 6. Effect of PO₄³⁻ on MaxNuclease endonuclease activity. PO₄³⁻ strongly inhibits the MaxNuclease endonuclease activity. The optimum PO₄³⁻ concentration is between 0 and 10mM.

Effect of Denaturant on Enzyme Activity




Figure 7. Effect of SDS on MaxNuclease endonuclease activity. SDS strongly inhibits the MaxNuclease endonuclease activity. 0.1% SDS inhibits nearly 90% of the activity.

Effect of Protein Precipitant Ammonium Sulfate on Enzyme Activity




Figure 8. Ammonium sulfate strongly inhibits MaxNuclease activity. Concentrations above 100 mM fully inhibit enzyme activity.

Effect of Surfactant (Tween 20) on Enzyme Activity




Figure 9. The effect of detergents on MaxNuclease endonuclease activity. The concentration of Tween 20 under 0.8% has no significant effect on MaxNuclease activity.

Temperature & Freeze/Thaw Stability




Figure 10. Temperature stability testing. MaxNuclease remains active when stored at 25°C/37°C for 7 days.




Figure 11. Freeze/thaw stability testing. Repeated freeze/thaws (up to 12) at -20°C / -80°C does not affect the performance.

High Purity: >99% via SEC-HPLC




Figure 12. The purity of MaxNuclease detected by SEC-HPLC is ≥ 99%.

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

