




PerfectStart[®] II Probe qPCR SuperMix UDG

Probe qPCR kit with high sensitivity and high amplification efficiency

Product Features

- › Three antibodies blocking technique, improves sensitivity, enhances specificity and amplification efficiency.
- › Specially optimized qPCR reaction buffer provides higher sensitivity and specificity.
- › dUTP/UDG used to avoid DNA carryover contamination and guarantee more accurate data.
- › Passive reference dyes are provided for different qPCR instruments to normalize well-to-well differences by pipetting errors or instrument limitations .
- › A wide range of applications, used for the detection of African swine fever virus, pseudorabies virus, Vibrio harveyi, Intestinal worm, and Vibrio parahaemolyticus.
- › Good stability which enables stable amplification when stored under repeated freezing and thawing, assembled reaction mix with primers and probes, room temperature, 37°C and other conditions.

Multiplex PCR Detection



Using Bovine Viral Diarrhea Virus (BVDV) cDNA, Porcine Blue-ear Disease Virus (PRRSV) cDNA, African Swine Fever Virus (ASFV) plasmid, and Pseudorabies virus (PRV) gDNA mixture as templates for quantitative detection.


The results show that TransGen product can perform 4-plex detection.

The best for life science


Distributed by:

CliniSciences Group

High Sensitivity and High Amplification Efficiency






qPCR was performed targeting different concentrations of PRRSV cDNA (10 pg, 0.1 pg, 0.01 pg) as input template, using products of TransGen, Company V, and Company T. The results show that the sensitivity of TransGen product can reach 0.1 pg, and the amplification effect is better than that of Company V and Company T.



qPCR was performed targeting different concentrations of Bovine Viral Diarrhea Virus (BVDV) cDNA (100 pg, 1 pg, 0.1 pg) as input template, using products of TransGen, Company V, and Company T. The results show that the sensitivity of TransGen product can reach 0.1 pg, and the amplification effect is better than that of Company V and Company T.



Wide Range of Applications



TransGen product can successfully detect different concentrations (10⁷ copies/ml-10⁴ copies/ml, 10-fold dilution) of the DNA of Vibrio harveyi, Enterocytozoon Hepatopenaei (EHP) and Vibrio parahaemolyticus.



Good Stability

Lot-to-lot Stability (Singleplex)




Singleplex qPCR was performed, using different lots of TransGen product. The results show that TransGen product can guarantee stable lot-to-lot repeatability.

Lot Stability (Quadruple Multiplex)




Quadruple multiplex qPCR was performed, using different lots of TransGen product. The results show that TransGen product can guarantee stable lot-to-lot repeatability.

Freeze-thaw Stability (Singleplex)




Singleplex qPCR was performed using TransGen product after 5, 10, 15 and 20 freeze-thaw cycles. The results show that the performance of TransGen product will not be affected after repeated freezing and thawing, and stable amplification can still be achieved.

Freeze-thaw Stability (Quadruple Multiplex)




Quadruple multiplex qPCR was performed using TransGen product after 5, 10, 15 and 20 freeze-thaw cycles. The results show that the performance of TransGen product will not be affected after repeated freezing and thawing, and stable amplification can still be achieved.

Assembled Reaction Mix Stability (Singleplex)




The primer and probes for a single detection gene were assembled with TransGen product, and qPCR detection was performed after 7 days of treatment at different temperatures. The results show that the TransGen assembled reaction mix can still be stably amplified after being stored at different temperatures.

Assembled Reaction Mix Stability (Double Multiplex)



The primers and probes for double detection genes were assembled with TransGen product, and qPCR detection was performed after 7 days of treatment at different temperatures. The results show that the TransGen assembled reaction mix can still be stably amplified after being stored at different temperatures.

Assembled Reaction Mix Stability (Triplex)



The primers and probes for triplex detection genes were assembled with TransGen product and qPCR detection was performed after 7 days of storage at different temperatures for 7 days. The results show that the TransGen assembled reaction mix can still be stably amplified after being stored at different temperatures.

Product Name	Cat. No.	Specification
PerfectStart® II Probe qPCR SuperMix UDG	AQ712-01	1 ml
	AQ712-02	5×1 ml
	AQ712-03	15×1 ml

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences Limited
Address: Marktgasse 18 8302 Kloten -
Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

