

Technical Specifications **ATR-001**

Principle :	Optical Imaging
Test Format :	Dipstick, Cassette or Panel
Measurement :	Qualitative/Semi-Quantitative/Quantitative
Read Time :	< 10 Seconds
Power Source :	5V DC 300mA
Dimensions(L*W*H):	155*119*111 mm
Weight :	440 g
OS :	Windows-based Proprietary Program
Connecting Ports:	USB 2.0
Analyte Information :	QR Code



Health Canada

EU IVDR Registered

CE Marked

Technical Specifications **ATR-200**

Principle:	Optical Imaging
Test Formats:	Cassette
Touch Screen:	5" LCD touch screen
Measurement:	Qualitative/Quantitative
Read Time:	< 60 seconds
Power Source:	AC 100-240 V / 50/60 Hz, DC 15 V
Dimensions (L*W*H):	270*127*115 mm
Weight:	< 1.2 kg
Memory:	8000 records
Connecting Ports:	RS232, USB
Printer:	Built in thermal printer



EU IVDR Registered

CE Marked

Technical Specifications **ATR-300**

Principle :	Optical Imaging and Image Processing
Test Format :	Dipstick, Cassette and Panel
Display:	7" LCD Touch Screen
Measurement :	Qualitative/Semi-Quantitative/Quantitative
Read Time :	5 Seconds
Power Source :	DC 12V, 2A
Dimensions(L*W*H) :	221 * 215 *165 mm
Weight :	< 1.5 kg
Data Storage :	Up to 3,000 Records
Connecting Ports :	LAN Port, 4 USB Ports
Printer :	Built-in thermal printer



EU IVDR Registered

CE Marked

Technical Specifications **ACR-100**

Principle :	Photoelectric Conversion
Test Format :	Cup
Display:	7" LCD Touch Screen
Measurement :	Qualitative
Read Time :	1 - 2.5 Minutes
Power Source :	DC 12V, 2A
Dimensions(L*W*H) :	221 * 215 *165 mm
Weight :	< 1.5 kg
Connecting Ports :	LAN Port, 4 USB Ports
Printer :	Built-in thermal printer



Health Canada

EU IVDR Registered

CE Marked

Hanazhou AllTest Biotech Co.,Ltd.

ALL TEST [®] Distributed by:
CliniSciences Group

LATERAL FLOW READER



Lateral Flow (LF) Readers
Rapid Tests For LF Readers
Cup Reader
DOA Rapid Tests For Cup Reader

Version:1 (English)

Quantitative Result Rapid Tests

PCT	Calprotectin	Ferritin	AMH	TSH
Pepsinogen	FOB	PSA	D-dimer	NGAL

Semi-Quantitative Result Rapid Tests

CRP	PCT	hCG	FSH	cTnI
cTnT	CK-MB	Myoglobin	D-dimer	FOB
Tetanus	TSH	Vitamin D		

Qualitative Result Rapid Tests

PCT	Syphilis	Cardiac Troponin I	Dengue IgG/IgM and NS1 Combo
CRP	Adeno Ag	Calprotectin	Rota and Adeno Combo
hCG	MONO	FOB	HBsAg/HCV/HIV/Syphilis Combo
LH	Dengue IgG/IgM	Legionella pneumophila	COVID-19 and Influenza A+B Ag Combo
fFN	Dengue NS1	Mycoplasma pneumoniae Ag	Noro/Rota/Adeno Combo
iGFBP-1	COVID-19 Ag	S. pneumoniae Ag	RSV/Flu A+B Combo
Chlamydia	COVID-19 IgG/IgM	Clostridium difficile GDH	S. pneumoniae and L. pneumophila Combo
H. pylori Ab	Strep A	D-dimer	Flu A+B+H1N1 Combo
H. pylori Ag	HbA1c	HBV Combo	
Flu A	IgE	Clostridium difficile GDH+Toxin A+Toxin B Combo	
Flu A+B	PSA	SARS-CoV-2/Influenza A+B/RSV Antigen Combo	
HBsAg	NT-proBNP	SARS-CoV-2/Flu A+B/RSV/Adeno Ag Combo	
HCV	Myoglobin	SARS-CoV-2/Flu A+B/RSV/Adeno/M.pneumoniae Ag Combo	
HIV 1.2.O	CK-MB	Adenovirus / RSV /Flu A+B Combo	

DOA Qualitative Result Rapid Tests


Oral Fluid			
Amphetamine (AMP)-50 ng/mL	Benzodiazepines (BZO)-30 ng/mL	Methadone (MTD)-30 ng/mL	
Barbiturate (BAR)-50 ng/mL	Cocaine (COC)-20 ng/mL	Opiates (OPI)-40 ng/mL	
Buprenorphine (BUP)-10 ng/mL	Methamphetamine (MET)-50 ng/mL	Phencyclidine (PCP)-10 ng/mL	
		Synthetic Marijuana (K2)-25 ng/mL	
Urine			
Alcohol (ALC)-0.02%	Ketamine (KET)-1000 ng/mL	Tricyclic Antidepressants (TCA)-1000 ng/mL	
Acetaminophen (ACE)-5000 ng/mL	Ecstasy (MDMA)-500 ng/mL	6-Monoacetylmorphine (6-MAM)-10 ng/mL	
Amphetamine (AMP)-1000 ng/mL	Tenamfetamine (MDA)-500 ng/mL	Lysergic Acid Diethylamide (LSD)-20 ng/mL	
Barbiturate (BAR)-300 ng/mL	Zolpidem (ZOL)-50 ng/mL	α-Pyrrolidinovalerophenone (α-PVP)-1000 ng/mL	
Buprenorphine (BUP)-10 ng/mL	Marijuana (THC)-50 ng/mL	Methylenedioxypyrovalerone (MDPV)-1000 ng/mL	
Benzodiazepines (BZO)-300 ng/mL	Tramadol (TML)-100 ng/mL	Methylphenidate (MPD)-300 ng/mL	
Clonazepam (CLO)-400 ng/mL	Opiates (OPI)-2000 ng/mL	Methamphetamine (MET)-1000 ng/mL	
Cocaine (COC)-300 ng/mL	Oxycodone (OXY)-100 ng/mL	Morphine (MOP)-300 ng/mL	
Cotinine (COT)-200 ng/mL	7-Aminoclonazepam (7-ACL)-300 ng/mL	Methaqualone (MQL)-300 ng/mL	
Diazepam (DIA)-300 ng/mL	Phencyclidine (PCP)-25 ng/mL	Methadone (MTD)-300 ng/mL	
Ethyl Glucuronide (ETG)-500 ng/mL	Propoxyphene (PPX)-300 ng/mL	Ethylenediamine-dimethylphosphinic acid (EDDP)-100 ng/mL	
Fentanyl (FYL)-20 ng/mL	Synthetic Marijuana (K2)-50 ng/mL		

Test Report for LF Reader

Test Report

1) Product
Analyte: Calprotectin Lot: CAL23020016E
Exp. Date: 2025-01

2) Patient
ID: Birthdate: 2000/1/1 Name: Gender: Male

3) Result

Item Calprotectin Result 339.5 ug/g (Positive)


4) Testing Site
Name: Address: Phone: Operator: admin
Test Time: 2023-03-23 15:15:29 Print Time: 2023-03-23 15:16:28

Quantitative Result

Test Report

1) Product
Analyte: D-dimer Lot: DDM22070006E
Exp. Date: 2024-06

2) Patient
ID: Birthdate: 2000/1/1 Name: Gender: Male

3) Result

Item D-dimer Result 2500-5000ng/mL (++)


4) Testing Site
Name: Address: Phone: Operator: admin
Test Time: 2023-03-23 15:15:14 Print Time: 2023-03-23 15:16:37

Semi-Quantitative Result

Test Report

1) Product
Analyte: H. pylori Antigen Lot: HPG23020051E
Exp. Date: 2025-01

2) Patient
ID: Birthdate: 2000/1/1 Name: Gender: Male

3) Result

Item H. pylori Antigen Result Positive

4) Testing Site
Name: Address: Phone: Operator: admin
Test Time: 2023-03-23 15:16:08 Print Time: 2023-03-23 15:16:16

Qualitative Result

Test Items for Cup Reader (Qualitative)

Multi-Drug Rapid Test 1-Step Cup With/Without Adulteration (Specimen:Urine)

AB-PINACA (ABP)-10 ng/mL
Acetaminophen (ACE)-5000 ng/mL
7-Aminoclonazepam (7-ACL)-300 ng/mL
Alprazolam(ALP)-100 ng/mL
Barbiturate (BAR)-200 ng/mL, 300 ng/mL
Buprenorphine (BUP)-5 ng/mL,10 ng/mL
Cathine (CAT)-150 ng/mL
Caffeine(CAF)-1000 ng/mL
Carisoprodol (CAR)-2000 ng/mL
Cannabinol (CNB)-500 ng/mL
Carfentanyl(CFYL)-250 ng/mL, 500 ng/mL
Citalopram(CIT)-500 ng/mL
Clonazepam (CLO)-150 ng/mL, 400 ng/mL
Cocaine (COC)-100 ng/mL, 150 ng/mL, 300 ng/ML
Diazepam (DIA)-300 ng/mL
Ethyl Glucuronide (ETG)-500 ng/mL, 1000 ng/mL
Fluoksetine (FKET)-1000 ng/mL
Fluoxetine (FLX)-500 ng/mL
Gabapentin(GAB)-2000 ng/mL
Kratom(KRA)-300 ng/mL
6-Monoacetylmorphine(6-MAM)-10 ng/mL
Methcathinone(MCAT)-500 ng/mL
Ecstasy (MDMA)-500 ng/mL, 1000 ng/mL
Tenamfetamine (MDA)-500 ng/mL
Mephedrone HCl (MEP)-100 ng/mL
Mescaline (MES)-100 ng/mL, 300 ng/mL
Morphine (MOP)-100 ng/mL, 300 ng/mL
Methylphenidate(MPD)-300 ng/mL, 1000 ng/mL
Meperidine (MPRD)-100 ng/mL
Methaqualone (MQL)-300 ng/mL
Methadone (MTD)-200 ng/mL, 300 ng/mL
Mirtazapine (MTZ)-500 ng/mL
Quetiapine(QTP)-1000 ng/mL

Methamphetamine (MET)-300 ng/mL, 500 ng/mL, 1000 ng/mL
N,N-Dimethyltryptamine (NND)-1000 ng/mL
Amphetamine (AMP)-300 ng/mL, 500 ng/mL, 1000 ng/mL
Benzodiazepines (BZO)-100 ng/mL, 200 ng/mL, 300 ng/mL
Cotinine (COT)-10 ng/mL, 50 ng/mL, 100 ng/mL, 200 ng/mL, 500 ng/mL
Ethylenediamine-dimethylphosphinic acid (EDDP)-300 ng/mL, 100 ng/mL
Scopolamine(SCOP)-500 ng/mL
Synthetic Marijuana (K2)-30 ng/mL, 50 ng/mL
Opiates (OPI)-2000 ng/mL,300 ng/mL
Oxycodone (OXY)-100 ng/mL
Olanzapine (OZP)-1000 ng/mL
Papaverine (PAP)-500 ng/mL
Phencyclidine (PCP)-25 ng/mL
Pregabalin (PGB)-500 ng/mL
Propoxyphene (PPX)-300 ng/mL
Risperidone (RPD)-150 ng/mL
UR-144 -25 ng/mL
Zaleplon(ZAL)-100 ng/mL
Zolpidem(ZOL)-50 ng/mL
Zopiclone(ZOP)-50 ng/mL
Tapentadol (TAP)-1000 ng/mL
Tilidine (TLD)- 50 ng/mL
Tramadol (TML)-100 ng/mL, 200 ng/mL, 300 ng/mL
Tropicamide(TRO)-350 ng/mL
Trazodone(TZD)-200 ng/mL
Tricyclic Antidepressants (TCA)-500 ng/mL,1000 ng/mL
Fentanyl (FYL)-10 ng/mL, 20 ng/mL, 100 ng/mL, 200 ng/mL
Ketamine (KET)-100 ng/mL, 300 ng/mL, 500 ng/mL, 1000 ng/mL
Lysergic Acid Diethylamide (LSD)-10 ng/mL,20 ng/mL, 50 ng/mL
Methylenedioxypyrovalerone (MDPV)-300 ng/mL, 500 ng/mL, 1000 ng/mL
α-Pyrrolidinovalerophenone (α-PVP)-300 ng/mL, 500 ng/mL, 1000 ng/mL
Marijuana (THC)-20 ng/mL, 25 ng/mL, 50 ng/mL, 150 ng/mL, 300 ng/mL

TEST REPORT

1) CUP
Ref.: DIA-R1137-A1
Lot: DOA23030092E
Exp.: 2025-02

2) PATIENT
ID:
Name:
Age: 0
Gender:

3) RESULT

ANALYTE	RESULT
CRE	NOR
PH	NOR
OPI	NEG
COC	NEG
MTD	NEG
TML	NEG
BAR	NEG
BZO	NEG
MET	NEG
PPX	NEG
AMP	NEG
KET	NEG
THC	NEG
K2	NEG
PGB	NEG

Tested: 2023-03-23 15:32:22
Printed: 2023-03-23 15:33:18

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkeudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

