

## ❖ Veterinary Urinalysis Reagent Dipsticks




**Specimen:** Urine  
**Measurement:** Qualitative / Semi-quantitative  
**Reading Time:** 30-120 Seconds  
**Shelf life:** 24 Months  
**Storage Temperature:** 2-30 °C  
**Analytes:** Ascorbic Acid (ASC), Glucose (GLU), Bilirubin (BIL), Ketone (KET), Blood (BLO), Leukocytes (LEU), Nitrite (NIT), pH, Specific Gravity (SG), Protein (PRO), Calcium (CA), Urobilinogen (URO), Creatinine (CRE), Microalbumin (ALB)

## ❖ Veterinary Urine Analyzer VUR-200




**Principle:** Photoelectric Colorimetry  
**Testing Speed :** 60 tests/hour  
**Data Memory:** 10,000 Results  
**Strip Incubation Time:** 1 minute  
**Light Wavelength:** 470nm, 525 nm and 625 nm  
**Display :** 3.5" LCD touch screen  
**Printer:** Built-in thermal printer  
**Connectivity:** RS232 port ,USB 2 ports  
**Working Condition:** 10-30°C, RHs 80% , 76-106kPa  
**Working Voltage:** AC100-240V, 50/60Hz, DC 5.85V±0.45V (AA×4)  
**Weight:** 600g  
**Dimensions(L/W/H):** 200\*137\*52 mm

## ❖ Veterinary Test Cup (Feces)



Cat. No.	Product	Specimen
VIFC-507	Feline Calicivirus Ag	Secretions
VIAIH-507	Avian Influenza Virus H5 (AIV H5) Ag	Laryngotracheal secretion/Feces
VIPV-607	Canine Parvovirus Ag	Feces/Vomit
VICV-607	CCV Ag	Feces/Vomit
VIFP-607	Feline Panleucopenia Virus Ag	Feces/Vomit
VIFCO-607	Feline Corona Virus Ag	Feces
VIRO-607	Rotavirus Ag	Feces
VIGA-607	Giardia Lambliia Ag	Feces
VITO-607	Toxoplasma Gondii Ag	Feces
VINDV-607	Newcastle Disease Virus Ag	Feces/Secretions
VIALV-607	Avian Leukosis Ag	Meconium/Cloaca Swab
VICR-607	Cryptosporidium Ag	Feces
VIDP-527	CDV+CPV Ag Combo	Secretions+Feces/Vomit
VIPC-627	CPV+CCV Ag Combo	Feces/Vomit
VIPG-627	CPV Ag+Giardia Lambliia Ag Combo	Feces
VIFPG-627	FPV+Giardia Lambliia Ag Combo	Feces
VIFF-627	FCoV Ag+FPV Ag	Feces
VIPCG-637	CPV + CCV+Giardia Ag Combo	Feces
VIFFG-637	FCoV+FPV+Giardia Lambliia Ag Combo	Feces
VIRCE-647	Rotavirus+Cryptosporidia+Coronavirus+E. Coli Ag Combo	Feces

(Pack: 5T/Format: Cup)



## ❖ Fluorescence Immunoassay Test

Cat. No.	Product	Specimen	Cut-off	Test Range	Result time	Sensitivity/Specificity/Accuracy	Sampling volume
VF-CPL-302-A	Canine Pancreatic Lipase (cPL) Ag	S/P	200 ng/mL	50 ng/mL-2000 ng/mL	15min	Sensitivity:>99.99%/Specificity:97.37%/Accuracy:98.68%	S/P: 10 µL
VF-FPL-302-A	Feline Pancreatic Lipase (fPL) Ag	S/P	3.6 ng/mL	2 ng/mL-50 ng/mL	15min	Sensitivity:98.18%/Specificity:98.00%/Accuracy:98.10%	S/P: 20 µL
VF-CNBNP-302-A	Canine NT-proBNP Ag	S/P	900pmol/L	500 pmol/L-10000 pmol/L	15min	Sensitivity:91.67%/Specificity:96.67%/Accuracy:94.17%	S/P: 100 µL
VF-CYSC-302-A	Canine Cystatin C (cCysC) Ag	S/P	1.2 mg/L	0.1 mg/L-10 mg/L	10min	Sensitivity:96.67%/Specificity:95.00%/Accuracy:96.00%	S/P: 10 µL
VF-COR-302-A	Cortisol (COR) Ag	S/P	/	1µg/dL-30µg/dL	15min	Sensitivity:94.83%/Specificity:94.23%/Accuracy:95.00%	S/P: 50 µL
VF-FYSC-302-A	Feline Cystatin C (fCys C) Ag	S/P	7 mg/dL	2 mg/dL-50 mg/dL	10min	Sensitivity:90.38%/Specificity:95.16%/Accuracy:92.98%	S/P: 20 µL
VF-CT4-302-A	Canine Thyroxine (cT4)	S/P	/	0.5 µg/dL-10 µg/dL	15min	Sensitivity:94.00%/Specificity:98.00%/Accuracy:96.00%	S/P: 50 µL
VF-FT4-302-A	Feline Thyroxine (fT4)	S/P	/	0.5 µg/dL-10 µg/dL	15min	Sensitivity:96.36%/Specificity:94.00%/Accuracy:95.24%	S/P: 50 µL
VF-CP4-302-A	Canine Progesterone (Cp4)	S/P	/	1 ng/mL-30 ng/mL	10min	Sensitivity:91.67%/Specificity:96.67%/Accuracy:94.17%	S/P: 75 µL
VF-CDD-302-A	Canine D-Dimer	S/P	0.5 mg/L	0.1 mg/L-10 mg/L	10min	Sensitivity:90.00%/Specificity:92.31%/Accuracy:91.18%	S/P: 20 µL
VF-CTSH-302-A	Canine Thyroid Stimulating Hormone (cTSH)	S/P	10 mg/L	0.25 ng/mL-5 ng/mL	15min	/	S/P: 75 µL
VF-CCR-402-A	Canine C-Reactive Protein (cCRP)	WB/S/P	0.5 mg/mL	5 mg/L-200 mg/L	5min	Sensitivity:95.00%/Specificity:92.50%/Accuracy:94.00%	WB: 6 µL; S/P: 3 µL
VF-FSA-402-A	Feline Serum Amyloid A (fSAA) Ag	WB/S/P	2mg/L	2 mg/L-200 mg/L	10min	Sensitivity:98.28%/Specificity:98.18%/Accuracy:98.23%	WB: 5 µL; S/P: 2.5 µL
VF-FLV-402-A	Feline Leukemia Virus (FeLV) Ag	WB/S/P	1.5 IU/mL	1.0 IU/mL-200 IU/mL	10min	Sensitivity:94.00%/Specificity:>99.99%/Accuracy:97.00%	WB: 10 µL; S/P: 5 µL
VF-FIV-402-A	Feline Immunodeficiency (FIV) Ab	WB/S/P	2.0 IU/mL	1.0 IU/mL-20 IU/mL	10min	Sensitivity:98.00%/Specificity:98.28%/Accuracy:98.15%	WB: 10 µL; S/P: 5 µL
VF-HW-402-A	Canine/Feline Heartworm (CHW/FHW) Ag	WB/S/P	800 IU/mL	30-20000 IU/mL	15min	Sensitivity:96.15%/Specificity:98.75%/Accuracy:98.11%	WB: 40 µL; S/P: 20 µL
VF-FCO-602-A	Feline Coronavirus (FCoV) Ag	Feces	10 IU/mL	5-1000 IU/mL	10min	Sensitivity:94.55%/Specificity:98.00%/Accuracy:96.19%	Moderate amount
VF-CPG-602-A	Canine Parvovirus (CPV) Ag	Feces	220 IU/mL	40-20000 IU/mL	10min	Sensitivity:97.14%/Specificity:97.50%/Accuracy:97.27%	Moderate amount
VF-FPG-602-A	Feline Panleucopenia Virus (FPV) Ag	Feces	100 IU/mL	50-20000 IU/mL	10min	Sensitivity:96.15%/Specificity:98.00%/Accuracy:97.06%	Moderate amount
VF-CCO-602-A	Canine Coronavirus (CCV) Ag	Feces	5 IU/mL	1-100 IU/mL	10min	/	Moderate amount
VF-RBV-802-A	Rabies Virus Ag	Saliva	100 IU/mL	30-100000 IU/mL	15min	/	Moderate amount
VF-CPDA-335-A	CPV/CDV/CAV Ab Combo	S/P	/	1 U- 6 U	15min	/	S/P: 30 µL
VF-FPHC-335-A	FPV/FHV/FCV Ab Combo	S/P	/	1 U- 6 U	15min	/	S/P: 30 µL

(Pack: 10T/Format: Cassette)

## ❖ Fluorescence Immunoassay Analyzer VFR-100

**Test Method :** Fluorescence Immunoassay  
**Test Format :** Cassette  
**Measurement :** Quantitative, Qualitative  
**Read Time :** < 20 Seconds  
**Display :** 7" LCD touch screen  
**Memory :** 10000 Records  
**Power Source:** 100-240VAC,50/60Hz  
**Dimensions (L\*W\*H):** 220\*190\*120 mm

**Weight:** 2 kg  
**Storage Condition:** 5-45°C, RH ≤85%  
**Working Condition:** 10-30°C, RH ≤85%  
**Printer:** Built-in thermal printer  
**Connectivity:** USB 2 ports, LAN port, RS232 COM port (WLAN is an optional function)  
**Scanner:** Laser barcode scanner with USB port  
**Compatibility:** LIS



## ❖ Lateral Flow (LF) Reader Plus VTR-200

**Principle:** Optical Imaging and Image Processing  
**Test Format:** Cassette  
**Touch Screen:** 7 Inches  
**Measurement:** Qualitative, Semi-quantitative, Quantitative  
**Read Time:** 5 Seconds  
**Power Source:** DC 12V, 2A  
**Dimensions (L\*W\*H):** 221 \* 215 \*165 mm  
**Weight:** < 1.5 kg  
**Data Storage:** Up to 3,000 Records  
**Connecting Ports:** LAN Port, 4 USB Ports,WIFI  
**Printer :** Built-in thermal printer



# VETERINARY TEST



Farm Animals (Live Stock)

Companion Animals (Pet)

Rapid Test

Dry Biochemistry

Fluorescence Immunoassay

Distributed by:

CliniSciences Group

Version:2 (English)





# CliniSciences Group

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oessterreich@clinisciences.com](mailto:oessterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [denmark@clinisciences.com](mailto:denmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Biotrend Chemikalien GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [info@biotrend.com](mailto:info@biotrend.com)  
Web: <https://www.biotrend.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Krajenhoffstraat 137A,  
1018RG Amsterdam, - Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences AS  
Address: c/o MerVerdi Munkerudtunet 10  
1164 Oslo - Norway  
Telephone: +47 21 988 882  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 916 750 700  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences AG  
Address: Fracht Ost Flughafen Kloten  
CH-8058 Zürich - Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: Biotrend Chemicals LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [info@biotrend-usa.com](mailto:info@biotrend-usa.com)  
Web: <https://www.biotrend-usa.com>

