

Product datasheet MON4044

MONOSAN[®]

Mouse anti-CD102, clone B-T1 (Monoclonal)

Clone no. B-T1

MONOSAN

Product name Mouse anti-CD102, clone B-T1 (Monoclonal)

Host Mouse

Applications WB,IF,FC,IP

Species reactivity human

Conjugate Purified

Immunogen ICAM-2 transfected CHO cells.

Isotype IgG1

Clonality Monoclonal

Clone number B-T1

Size 0.2 mg

Concentration 1.0 mg/ml

Format -

Storage buffer PBS

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD102, clone B-T1 (Monoclonal)

Clone no. B-T1

MONOSAN

Additional info

Mouse anti Human CD102 antibody, clone B-T1 recognizes human Intercellular adhesion molecule 2, also known as CD102 or ICAM-2. CD102 is a 275 amino acid ~55-65 kDa single pass type-1 transmembrane glycoprotein containing two Ig-like C2-type domains. Mouse anti Human CD102 antibody, clone B-T1 inhibits cell adhesion (Xie et al. 1995) and T cell activation and also recognizes soluble ICAM-2.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES