



Normal Sheep Serum

CatNo.: ENS9010

Lot Number: See Label

Size: CatNo.: ENS9010-2; 2ml
CatNo.: ENS9010-10; 10ml

Application: Immunhistology – Frozen
Immunhistology – Paraffin
Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures.

Suggested Working Dilution For blocking purpose, use 2.5 – 5% (v/v) solution.

General notes: Normal sheep serum can be used for blocking non-specific conserved sequence and/or Fc receptor binding.
Typically, the serum used for blocking is from a different species than the species in which the primary antibody was raised. Often the blocking serum is from the species in which the secondary antibody was raised.

Collection: When searchers choose their serum, an important factor that should be taken into consideration is the source, which also emphasizes the traceability of the serum. Our system of vertical integration allows us to be certain of the origin and traceability of our serum. The serum is imported and treated in agreement with the European regulations. It is heat treated at 56°C for 30 minutes.

Protein Concentration: 56 mg/ml

Preservative: 0.08% NaN₃

Shipping: Ship at room temperature.

Storage: Store at 2-8°C for 12 month or at -20°C for longer.

Shelf Life: 12 months from date of dispatch.

Regulatory: For **in vitro research** use only, not for diagnostic or therapeutic use.

Health and Safety Information:

This product contains sodium azide (NaN₃): a substance which should be handled by trained staff only.

Manufacturer

BIOZOL
FIT FOR SCIENCE

BIOZOL

Diagnostica Vertrieb GmbH
Leipziger Straße 4
85386 Eching

Phone +49 (89) 3799 666 6
Fax +49 (89) 3799 666 99
E-Mail info@biozol.de
www biozol.de

