




Anti Mouse Dab 1 Polyclonal Antibody

Package Size	25 μ g (100 μ L / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20°C Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with His tagged protein containing PTB domain and tyrosine-phosphorylated domain by antigen affinity chromatography.
Working dilution	For Immunoprecipitation ;5 μ g/mL

Western blotting (After Immunoprecipitation)



Sample : Flag-Dab1 transfected cell lysate

- ① Immunoprecipitated with Anti Flag Antibody
- ② Immunoprecipitated with Anti Dab1 Antibody

Preparation of antibodies and instruction :

Mikoshiha K, Morimura T
Department of Basic Medical Sciences,
The Institute of Medical Science, The University of Tokyo



Anti Mouse Dab 1 Polyclonal Antibody

【Reference】

1. Gleeson JG. et al. : Trends Neurosci. 2000 Aug;23(8):352-9
2. Howell BW. et al. : EMBO. 1997 Jan 2;16(1):121-32
3. Ogawa M. et al. : Neuron. 1995 May;14(5):899-912
4. Sheldon M. et al. : Nature. 1997 Oct 16;389(6652):730-3
5. Utsunomiya-Tate N. et al. : Proc Natl Acad Sci USA. 2000 Aug 15;97(17):9729-34

Manufacturer



Medicinal Chemistry Pharmaceutical Co., Ltd.

Kobe Research Institute

7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-945-7075 FAX:+81-78-306-0694

URL:<https://soyaku.co.jp/english/> tech-kobe@soyaku.co.jp

Previous manufacturer



Trans Genic Inc.