Product datasheet MON4068



Mouse anti-CD89, clone MIP8a (Monoclonal)

Clone no. MIP8a MONOSAN

Product name Mouse anti-CD89, clone MIP8a (Monoclonal)

Host Mouse

Applications IP,FC,ELISA

Species reactivity human

Conjugate Purified

Immunogen Recombinant soluble human Fc alpha R.

lsotype lgG1

Clonality Monoclonal

Clone number MIP8a

Size 0.1 mg

Concentration 1.0 mg/ml

Format -

Storage buffer PBS with 0.09% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON4068



Mouse anti-CD89, clone MIP8a (Monoclonal)

Clone no. MIP8a MONOSAN

Additional info

Mouse anti Human CD89 antibody, clone MIP8a recognizes the human CD89 cell surface antigen, a 50-75 kDa cell surface glycoprotein that is also known as the IgA receptor (Fc alpha R).CD89 is expressed by peripheral blood monocytes and neutrophils.MIP8a blocks binding of IgA to the Fc alpha R, and also inhibits neutrophil phagocytosis of IgA complexes.

References

1. -

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES