

A19220

Leader in Biomolecular Solutions for Life Science



Rabbit Anti-human immunoglobulin G mAb

Catalog No.: A19220

Basic Information

Observed MW

Refer to figures

Calculated MW

150kDa

Category

SMab Recombinant Monoclonal
Antibody

Applications

ELISA, Specific methodologies have not
been tested.

Cross-Reactivity

Human

Background

The Rabbit Anti-human immunoglobulin M mAb is developed from the heavy and light chain of all human immunoglobulin G.

Recommended Dilutions

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Purified Human IgG(H&L chain)

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

>90% pure (SDS-PAGE)

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: 10mM Sodium Phosphate, 150 mM Sodium Chloride, pH 7.4