

## Pipette Controller Cat# NB-12-6005-UK-PLUG



Technical Data	
Power	2 x 1.2 V / 1600 mAh NiMH batteries
Battery Service Life	Eight hours continuous use
Battery Charging Time	2 - 3 hours to charge fully
Volumetric Range	Plastic or Glass pipette from 0.1 ml up to 100 ml
Safety System	Replaceable hydrophobic 0.45 micron PTFE filter, Safety check valve
Standard Accessories	2 NiMH Batteries, Battery charger, Desktop Stand, Wall mounting stand, Hydrophobic PTFE membrane filter 0.45 micron 2Pcs, Hydrophobic PTFE membrane filter 0.20 micron 1Pc