

## Agarose Electrophoresis Grade

Cat # NB-42-00145-100g    Size: 100g

Cat # NB-42-00145-500g    Size: 500g

Cat # NB-42-00145-1Kg    Size: 1kg

### Product Information

Gel Strength	1.5% >100g/cm <sup>2</sup> - 2200g/cm <sup>2</sup>
Gelling temp	36 ± 1.5C
Melting temp	88 ± 1.5C
Ash Content	<0.45%
Sulphate	<0.1%
Clarity	<4 NTU
DNA/RNase	Negative.
Use Rates	0.75% - 2%.
EEO	0.14 - 0.16
CAS	9012-36-6
MDL Number	MFCD00081294
EINECS Number	232-731-8

### Storage

Store at room temperature.

*For reference only  
For Research Use Only. Not for Diagnostic or Therapeutic Use.*