

Product datasheet PS110P

MONOSAN[®]

Polyclonal Rabbit anti-Rat IgG, IgM and IgA, HRP

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Rat IgG, IgM and IgA, HRP
Host	Rabbit
Applications	WB, ELISA (1:1000), IHC
Species reactivity	rat
Conjugate	HRP
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS110P

MONOSAN[®]

Polyclonal Rabbit anti-Rat IgG, IgM and IgA, HRP

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/ml in the buffer (0.05M sodium phosphate pH 7.5, 0.1 M NaCl) containing 50% glycerol and 5 mg/ml bovine serum albumin. Rat immunoglobulins IgG, IgM, IgA

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES
