

A18859

Leader in Biomolecular Solutions for Life Science



IAA3 Rabbit pAb

Catalog No.: A18859

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Arabidopsis thaliana

Background

SHY2/IAA3 regulates multiple auxin responses in roots. It is induced rapidly by IAA, and has been shown to be phosphorylated by oat phytochrome A in vitro.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

839570

Swiss Prot

Immunogen

Recombinant protein of Arabidopsis thaliana IAA3.

Synonyms

F19P19.32; F19P19_32; IAA3; indole-3-acetic acid inducible 3; SHORT HYPOCOTYL 2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.