ABclonal®

Recombinant Mouse IL-22 Protein

Catalog No.: RP01617 Recombinant

Sequence Information

Species Gene ID Swiss Prot HEK293 cells 50929 Q9JJY9

Tags C-hFC

Synonyms

Interleukin-22; IL-22; Cytokine Zcyto18; IL-TIF; IL-D110; TIFa; TIFIL-23; ZCYTO18;IL22

Product Information

Source HEK293 cells Purification > 95% by SDS-

PAGE.

Endotoxin

<0.1EU/µg

Formulation

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

The cytokine interleukin-22 (IL-22) is a critical regulator of epithelial homeostasis. It has been implicated in multiple aspects of epithelial barrier function, including regulation of epithelial cell growth and permeability, production of mucus and antimicrobial proteins (AMPs), and complement production.

Basic Information

Description

Recombinant Mouse IL-22 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Leu34-Val179) of mouse IL-22 (Accession #NP_058667.1) fused with and a hFc tag at the C-terminus.

Bio-Activity


Measured by its ability to induce IL-10 secretion in COLO 205 human colorectal adenocarcinoma cells. The ED₅₀ for this effect is 150-600 ng/mL.

Storage


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.
 -80°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Mouse IL-22 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 60-68 kDa.



Recombinant Mouse IL-22 induce IL-10 secretion in COLO 205 human colorectal adenocarcinoma cells. The ED₅₀ for this effect is 150-600 ng/mL.