




Cloud-Clone Corp.

Distributed by:

CliniSciences Group




Cloud-Clone Antibodies

**You may not be aware of
four advantages!**

CLOUD-CLONE CORP. (CCC, USA), a research and experimental development in molecular biology and immunology oriented company headquartered in Houston, is the owner of the trade mark "Cloud-Clone".



1. Original Manufacturer for Antibodies



2. WB+IHC+IF Applications Featured Products


Based on more than ten years of research and development technology accumulation, Cloud-Clone has developed WB+IHC+IF applications featured products, which can be widely used in human, mouse, rat, rabbit, porcine, bovine, caprine and other species to assist scientific researchers in experimental research.



Fig. (A) WB analysis of extracts of tissues and cells of different species using TUBb antibody; (B)(C)(D) IHC analysis of different tissues using TUBb antibody; (E)(F) IF analysis of different cells using TUBb antibody

3. Antibodies can be used to prepare high-quality ELISA kits

In March 2020, CiteAb, the world-renowned biological reagent search engine, awarded Cloud-Clone as the “ELISA kit supplier to watch in 2020”. This is a great inspiration of Cloud-Clone. The number of citations of Cloud-Clone ELISA kits come out top in terms of various disease and cancer markers, cytokines and signaling pathways, etc.. In recent years, Cloud-Clone have been fully realized self-research, self-production and self-marketing. Antibodies and antibody pairs are the raw materials for ELISA kits, so ensuring the quality of antibodies is the most significant step.



4. Excellent literatures of Cloud-Clone Antibodies



Product No.	Product Name	Magazine	PubMed ID/doi
PAB030Mu01	Polyclonal Antibody to Protein Tyrosine Phosphatase Receptor Type C (CD45)	Cell	29958103
PAA102Hu01	Polyclonal Antibody to Matrix Metalloproteinase 7 (MMP7)	Science Translational Medicine	29167395
PAB755Mu01	Polyclonal Antibody to Interferon Regulatory Factor 4 (IRF4)	Blood	30002144
PAB153Hu01	Polyclonal Antibody to Leukocyte Immunoglobulin Like Receptor Subfamily B, Member 2 (LILRB2)	Nature Communications	29467366
PAA816Mu01	Polyclonal Antibody to Alpha-1-Acid Glycoprotein (a1AGP)	EMBO Molecular Medicine	28927749

Top 100 antibodies sales list:

Monoclonal antibodies	
Product No.	Product Name
MAA541Mu21	Monoclonal Antibody to Brain Natriuretic Peptide (BNP)
MAA077Hu22	Monoclonal Antibody to Interleukin 4 (IL4)
MAB120Hu22	Monoclonal Antibody to Arginase (ARG)
MAA182Ge21	Monoclonal Antibody to Hyaluronic Acid (HA)
MAA855Ra21	Monoclonal Antibody to Receptor Activator Of Nuclear Factor Kappa B Ligand (RANKL)
MAA927Ge21	Monoclonal Antibody to Histamine (HA)
MAB526Ge21	Monoclonal Antibody to Lipopolysaccharide (LPS)
MAB459Ge21	Monoclonal Antibody to Ovalbumin (OVA)
MAA071Hu22	Monoclonal Antibody to Interleukin 1 Alpha (IL1a)
MAC130Hu22	Monoclonal Antibody to Noggin (NOG)
MAA665Hu21	Monoclonal Antibody to Cross Linked C-Telopeptide Of Type I Collagen (CTXI)
MAL621Mu21	Monoclonal Antibody to Interleukin 1 Family, Member 9 (IL1F9)
MAB780Hu22	Monoclonal Antibody to Orosomucoid 2 (ORM2)
MAA139Hu22	Monoclonal Antibody to Tumor Necrosis Factor Related Apoptosis Inducing Ligand (TRAIL)
MAB758Hu22	Monoclonal Antibody to Ribonuclease A3 (RNASE3)
MAA857Hu21	Monoclonal Antibody to Clara Cell Protein 16 (CC16)
MAA957Hu21	Monoclonal Antibody to Procollagen I N-Terminal Propeptide (PINP)
MAG531Hu21	Monoclonal Antibody to Lysophosphatidylcholine Acyltransferase 3 (LPCAT3)
MAA573Hu21	Monoclonal Antibody to Procollagen III N-Terminal Propeptide (PIIINP)
MAA485Mi21	Monoclonal Antibody to N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP)
MAC000Hu22	Monoclonal Antibody to Haptoglobin Related Protein (HPR)
MAC489Hu21	Monoclonal Antibody to Ferroportin (FPN)
MAB555Hu22	Monoclonal Antibody to A Disintegrin And Metalloprotease 17 (ADAM17)
MAB980Hu24	Monoclonal Antibody to Interleukin 33 (IL33)
MAA819Ra21	Monoclonal Antibody to Parathyroid Hormone Related Protein (PTHrP)
MAA124Hu22	Monoclonal Antibody to Transforming Growth Factor Beta 1 (TGFb1)
MAA207Hu22	Monoclonal Antibody to Alanine Aminotransferase (ALT)
MAA462Ge21	Monoclonal Antibody to Cortisol (Cor)
MAA143Hu21	Monoclonal Antibody to Vascular Endothelial Growth Factor A (VEGFA)
MAA222Hu22	Monoclonal Antibody to Interferon Beta (IFNb)
MAB215Hu22	Monoclonal Antibody to Fibrinogen Beta (FGb)
MAC109Hu21	Monoclonal Antibody to Growth Differentiation Factor 3 (GDF3)
MAF551Hu21	Monoclonal Antibody to Lysyl Oxidase Like Protein 1 (LOXL1)

Polyclonal antibodies	
Product No.	Product Name
PAA788Hu01	Polyclonal Antibody to Programmed Cell Death Protein 1 Ligand 1 (PDL1)
PAG854Ge01	Polyclonal Antibody to Phosphatidylinositol Phosphate (PIP)
PAC614Ra01	Polyclonal Antibody to Macrophage Receptor With Collagenous Structure (MARCO)
PAN056Mu01	Polyclonal Antibody to T-Cell Immunoreceptor With Ig And ITIM Domains Protein (TIGIT)
PAB701Ge01	Polyclonal Antibody to Cholesterol (CH)
PAB366Ra01	Polyclonal Antibody to Cadherin 5 (CDH5)
PAG531Hu01	Polyclonal Antibody to Lysophosphatidylcholine Acyltransferase 3 (LPCAT3)
PAB934Mu01	Polyclonal Antibody to Leucine Rich Alpha-2-Glycoprotein 1 (LRG1)
PAN056Hu01	Polyclonal Antibody to T-Cell Immunoreceptor With Ig And ITIM Domains Protein (TIGIT)
PAA593Ra03	Polyclonal Antibody to Fibrillin 1 (FBN1)
PAB526Ge01	Polyclonal Antibody to Lipopolysaccharide (LPS)
PAS099Ge01	Polyclonal Antibody to Sialic Acid (SA)
PAB933Bo01	Polyclonal Antibody to Cathepsin S (CTSS)
PAB345Hu01	Polyclonal Antibody to Tetraspanin 30 (TSPAN30)
PAB197Hu01	Polyclonal Antibody to Cartilage Oligomeric Matrix Protein (COMP)
PAA182Ge01	Polyclonal Antibody to Hyaluronic Acid (HA)
PAJ089Mu01	Polyclonal Antibody to FOS Like Antigen 1 (FOSL1)
PAG515Mu01	Polyclonal Antibody to Chymase 1, Mast Cell (CMA1)
PAB305Mu01	Polyclonal Antibody to Complement Component 4b (C4b)
PAC994Hu01	Polyclonal Antibody to Glutathione Peroxidase 4 (GPX4)
PAB558Hu01	Polyclonal Antibody to Cytochrome P450 21B (CYP21B)
PAC141Hu01	Polyclonal Antibody to Collagen Type IV Alpha 5 (COL4a5)
PAB224Hu01	Polyclonal Antibody to Gastrin (GT)
PAE248Hu01	Polyclonal Antibody to Cysteine Dioxygenase I (CDO1)
PAB575Hu01	Polyclonal Antibody to Secondary Lymphoid Tissue Chemokine (SLC)
PAE914Ra01	Polyclonal Antibody to Sirtuin 4 (SIRT4)
PAA350Hu02	Polyclonal Antibody to Collagen Type I Alpha 1 (COL1a1)
PAF417Hu01	Polyclonal Antibody to Fatty Acid Desaturase 2 (FADS2)
PAH750Hu01	Polyclonal Antibody to Liver Expressed Antimicrobial Peptide 2 (LEAP2)
PAF553Mu01	Polyclonal Antibody to Lysyl Oxidase Like Protein 3 (LOXL3)
PAN021Hu01	Polyclonal Antibody to RNA Binding Motif Protein 20 (RBM20)
PAA817Ra03	Polyclonal Antibody to Haptoglobin (Hpt)

Distributed by:

CliniSciences Group

MAC022Ge21	Monoclonal Antibody to Reverse Triiodothyronine (rT3)
MAB214Hu21	Monoclonal Antibody to Aspartate Aminotransferase (AST)
MAC141Hu21	Monoclonal Antibody to Collagen Type IV Alpha 5 (COL4a5)
MAB758Hu21	Monoclonal Antibody to Ribonuclease A3 (RNASE3)
MAB526Ge22	Monoclonal Antibody to Lipopolysaccharide (LPS)
MAA156Hu21	Monoclonal Antibody to Carbohydrate Antigen 19-9 (CA19-9)
MAB353Ge21	Monoclonal Antibody to Advanced Glycation End Product (AGE)
MAC489Hu26	Monoclonal Antibody to Ferroportin (FPN)
MAC823Hu22	Monoclonal Antibody to Thymidine Kinase 1, Soluble (TK1)
MAA916Ge21	Monoclonal Antibody to Vitamin B6 (VB6)
MAA235Hu21	Monoclonal Antibody to Complement 1 Inhibitor (C1INH)
MAA056Ra21	Monoclonal Antibody to Interleukin 10 (IL10)
MAB337Hu21	Monoclonal Antibody to Elastin (ELN)
MAC128Hu22	Monoclonal Antibody to Gremlin 1 (GREM1)
MAA693Hu22	Monoclonal Antibody to Heat Shock 27kDa Protein 1 (HSPB1)
MAA859Hu21	Monoclonal Antibody to Mannose Associated Serine Protease 2 (MASP2)
MAH756Hu21	Monoclonal Antibody to Klotho Beta (KLB)
MAA055Hu22	Monoclonal Antibody to Insulin Like Growth Factor Binding Protein 4 (IGFBP4)
MAA079Rb21	Monoclonal Antibody to Interleukin 6 (IL6)
MAA748Ra21	Monoclonal Antibody to Trefoil Factor 2 (TFF2)
MAE676Hu26	Monoclonal Antibody to Regenerating Islet Derived Protein 3 Gamma (REG3g)
MAA428Hu21	Monoclonal Antibody to Neurokinin A (NKA)
MAA461Ge21	Monoclonal Antibody to Estradiol (E2)
MAA573Hu22	Monoclonal Antibody to Procollagen III N-Terminal Propeptide (PIIINP)
MAA991Mi28	Monoclonal Antibody to Ghrelin (GHRL)
MAB979Hu26	Monoclonal Antibody to Hepcidin (Hepe)
MAA133Ra21	Monoclonal Antibody to Tumor Necrosis Factor Alpha (TNFa)
MAA553Hu22	Monoclonal Antibody to Matrix Metalloproteinase 9 (MMP9)
MAA833Rb21	Monoclonal Antibody to Von Willebrand Factor (vWF)
MAS085Ra21	Monoclonal Antibody to Angiotensin 1-7 (Ang1-7)
MAA545Po21	Monoclonal Antibody to Immunoglobulin E (IgE)
MAA519Rb21	Monoclonal Antibody to Apolipoprotein A1 (APOA1)
MAB997Mu21	Monoclonal Antibody to Apolipoprotein A5 (APOA5)
MAB980Ra21	Monoclonal Antibody to Interleukin 33 (IL33)
MAA003Ge01	Monoclonal Antibody to Cyclic Adenosine Monophosphate (cAMP)
MAA788Hu22	Monoclonal Antibody to Programmed Cell Death Protein 1 Ligand 1 (PDL1)

PAA018Mu01	Polyclonal Antibody to Chemokine C-Motif Receptor 1 (XCR1)
PAH197Hu01	Polyclonal Antibody to Suprabasin (SBSN)
PAA908Ge01	Polyclonal Antibody to Melatonin (MT)
PAJ557Hu01	Polyclonal Antibody to Beclin 1 (BECN1)
PAA621Ge01	Polyclonal Antibody to Hydroxyproline (Hyp)
PAB154Mu02	Polyclonal Antibody to Fibrinogen Alpha (FGa)
PAH756Hu01	Polyclonal Antibody to Klotho Beta (KLB)
PAL621Mu01	Polyclonal Antibody to Interleukin 1 Family, Member 9 (IL1F9)
PAA642Hu01	Polyclonal Antibody to Aldolase A, Fructose Biphosphate (ALDOA)
PAB945Hu01	Polyclonal Antibody to Epiregulin (EREG)
PAD059Hu01	Polyclonal Antibody to Glucagon Like Peptide 2 (GLP2)
PAB550Hu01	Polyclonal Antibody to Poliovirus Receptor (PVR)
PAB978Mu01	Polyclonal Antibody to Complement Factor I (CFI)
PAA132Ra01	Polyclonal Antibody to Adiponectin Receptor 2 (ADIPOR2)
PAA843Mu01	Polyclonal Antibody to Gonadotropin Releasing Hormone (GnRH)
PAB558Mu01	Polyclonal Antibody to Cytochrome P450 21B (CYP21B)
PAB697Mu01	Polyclonal Antibody to Alpha-1-Antitrypsin (a1AT)
PAC914Hu01	Polyclonal Antibody to Fibroblast Growth Factor 12 (FGF12)
PAS089Ge01	Polyclonal Antibody to Deoxycholate (DC)
PAA751Hu01	Polyclonal Antibody to Programmed Cell Death Protein 1 (PDCD1)
PAB203Hu01	Polyclonal Antibody to Neurotensin (NT)
PAC611Hu01	Polyclonal Antibody to Zinc Finger Protein 10 (ZNF10)
PAB019Hu01	Polyclonal Antibody to Angiopoietin Like Protein 4 (ANGPTL4)
PAD696Mu01	Polyclonal Antibody to Diazepam Binding Inhibitor (DBI)
PAB406Mu02	Polyclonal Antibody to Lipopolysaccharide Binding Protein (LBP)
PAB611Hu02	Polyclonal Antibody to Neuraminidase (NEU)
PAA468Ra01	Polyclonal Antibody to Decay Accelerating Factor (DAF)
PAA980Ra01	Polyclonal Antibody to Endoglin (ENG)
PAB366Hu02	Polyclonal Antibody to Cadherin 5 (CDH5)
PAL060Hu01	Polyclonal Antibody to Parkinson Disease Protein 2 (PARK2)
PAB796Ra01	Polyclonal Antibody to Signal Transducer And Activator Of Transcription 2 (STAT2)
PAC038Ge01	Polyclonal Antibody to Diacylglycerol (DAG)
PAB573Mu01	Polyclonal Antibody to Mucosae Associated Epithelia Chemokine (MEC)
PAA252Hu01	Polyclonal Antibody to Apolipoprotein C1 (APOC1)
PAK522Ge01	Polyclonal Antibody to Bilirubin (Bb)
PAJ331Bo01	Polyclonal Antibody to Casein Kappa (CSN3)

MAA833Ra21	Monoclonal Antibody to Von Willebrand Factor (vWF)
MAB980Hu22	Monoclonal Antibody to Interleukin 33 (IL33)
MAB653Hu21	Monoclonal Antibody to Myostatin (MSTN)
MAB792Hu21	Monoclonal Antibody to S100 Calcium Binding Protein A8 (S100A8)
MAB793Hu22	Monoclonal Antibody to S100 Calcium Binding Protein A9 (S100A9)
MAA133Hu22	Monoclonal Antibody to Tumor Necrosis Factor Alpha (TNFa)
MAB166Hu22	Monoclonal Antibody to VGF Nerve Growth Factor Inducible (VGF)
MAA049Hu22	Monoclonal Antibody to Interferon Gamma (IFNg)
MAA087Hu22	Monoclonal Antibody to Monocyte Chemotactic Protein 1 (MCP1)
MAE787Hu22	Monoclonal Antibody to Annexin A1 (ANXA1)
MAA079Ra21	Monoclonal Antibody to Interleukin 6 (IL6)
MAA524Ra21	Monoclonal Antibody to Tissue Factor (TF)
MAA633Hu21	Monoclonal Antibody to Platelet Derived Growth Factor BB (PDGFBB)
MAA091Hu21	Monoclonal Antibody to Macrophage Derived Chemokine (MDC)
MAA835Mi21	Monoclonal Antibody to Corticotropin Releasing Hormone (CRH)
MAG962Hu21	Monoclonal Antibody to Tumor Protein, Translationally Controlled 1 (TPT1)
MAS085Mi21	Monoclonal Antibody to Angiotensin 1-7 (Ang1-7)
MAB877Hu22	Monoclonal Antibody to Forkhead Box Protein P3 (FOXP3)
MAA478Hu22	Monoclonal Antibody to Troponin I Type 3, Cardiac (TNNI3)
MAA726Hu22	Monoclonal Antibody to Transthyretin (TTR)
MAA920Ge21	Monoclonal Antibody to Vitamin D3 (VD3)
MAC174Hu21	Monoclonal Antibody to Neuregulin 4 (NRG4)
MAB995Hu22	Monoclonal Antibody to Hemojuvelin (HJV)
MAA049Gu21	Monoclonal Antibody to Interferon Gamma (IFNg)
MAB769Hu22	Monoclonal Antibody to S100 Calcium Binding Protein A6 (S100A6)
MAA563Hu21	Monoclonal Antibody to Interleukin 1 Beta (IL1b)
MAA757Hu22	Monoclonal Antibody to Epidermal Growth Factor Receptor (EGFR)
MAA816Hu22	Monoclonal Antibody to Alpha-1-Acid Glycoprotein (a1AGP)
MAB320Hu22	Monoclonal Antibody to Thymic Stromal Lymphopoietin (TSLP)
MAC147Hu21	Monoclonal Antibody to Collagen Type VII (COL7)
MAB547Hu22	Monoclonal Antibody to Indoleamine-2,3-Dioxygenase (IDO)

PAA253Hu01	Polyclonal Antibody to Maltase Glucoamylase, Intestinal (MGA)
PAB744Hu01	Polyclonal Antibody to Ubiquitin D (UBD)
PAC802Hu01	Polyclonal Antibody to Tensin 1 (TNS1)
PAH548Mu01	Polyclonal Antibody to Nexilin (NEXN)
PAJ226Hu02	Polyclonal Antibody to Biglycan (BGN)
PAA171Hu01	Polyclonal Antibody to Meprin A Alpha (MEP1a)
PAB468Mu01	Polyclonal Antibody to Myeloproliferative Leukemia Virus Oncogene (MPL)
PAC011Ra03	Polyclonal Antibody to Complement Factor B (CFB)
PAB503Hu01	Polyclonal Antibody to Dystrophin (DMD)
PAF959Hu01	Polyclonal Antibody to Utrophin (UTRN)
PAA038Mu01	Polyclonal Antibody to FMS Like Tyrosine Kinase 3 Ligand (Flt3L)
PAG905Hu01	Polyclonal Antibody to Vasorin (VASN)
PAA817Mu01	Polyclonal Antibody to Haptoglobin (Hpt)
PAC099Hu01	Polyclonal Antibody to Integrin Beta 6 (ITGb6)
PAC564Hu01	Polyclonal Antibody to Legumain (LGMN)
PAC598Hu01	Polyclonal Antibody to Vanin 1 (VNN1)
PAG608Hu01	Polyclonal Antibody to Galactosylceramidase (GALC)
PAB967Mu01	Polyclonal Antibody to Apolipoprotein A4 (APOA4)
PAD198Hu01	Polyclonal Antibody to Sequestosome 1 (SQSTM1)
PAB508Hu01	Polyclonal Antibody to Killer Cell Lectin Like Receptor Subfamily C, Member 2 (KLRC2)
PAB340Mi01	Polyclonal Antibody to Actin Beta (ACTb)
PAA934Hu01	Polyclonal Antibody to Peroxisome Proliferator Activated Receptor Alpha (PPARA)
PAD069Mu01	Polyclonal Antibody to Discs, Large Homolog 5 (DLG5)
PAG628Ra01	Polyclonal Antibody to Triggering Receptor Expressed On Myeloid Cells 2 (TREM2)
PAH755Ra01	Polyclonal Antibody to Kynurenine-3-Monooxygenase (KMO)
PAA866Bo01	Polyclonal Antibody to Parathyroid Hormone (PTH)
PAJ076Hu01	Polyclonal Antibody to Galactose-3-O-Sulfotransferase 1 (GAL3ST1)
PAC822Hu01	Polyclonal Antibody to Thromboxane Receptor (TP)
PAG152Hu01	Polyclonal Antibody to Small Ubiquitin Related Modifier Protein 2 (SUMO2)
PAB732Mu01	Polyclonal Antibody to C-Type Lectin Domain Family 10, Member A (CLEC10A)
PAD855Mu01	Polyclonal Antibody to Alanine Glyoxylate Aminotransferase 2 (AGXT2)
PAJ829Hu02	Polyclonal Antibody to Phosphoinositide-3-Kinase Catalytic Beta Polypeptide (PIK3Cb)

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: danmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mausier-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkeudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

