




We are one of the UK's leading suppliers of life science research solutions.

An outstanding, independent UK company offering quality products and first-class customer service to leading pharma, biotech, contract research organisations and universities worldwide.

We offer our 30 years of technical expertise and experience to provide traceable, customisable solutions through the three divisions of LSG: Life Science Production, Antibody Production Services and Diagnostic Kit Services.

All requirements are backed by our ISO 13485, ISO 9001:2015 and ISIA Traceability accreditations.



Technical Expertise and Industry Experience

With over 30 years of experience in antibody production, cell culture, diagnostics industries, and supporting the biopharmaceutical sector, we offer unparalleled technical expertise and industry insight. Our team is equipped to provide tailored guidance and support, covering project management, regulatory compliance, and more.

Complete Confidentiality

We are qualified with NDAs to ensure the confidentiality of your sensitive information.

Traceability

LSG is a member of the International Serum Industry Association (ISIA) and is ISIA Traceability Certified.

Quality Commitment

We are ISO 13485 and ISO 9001:2015 Quality Management System Certified.

Our Services



Custom Contract Services




Custom Media Formulation



Product Sourcing





Biospecimens



Bespoke Labelling



Antibody Production Services



Diagnostic Kit Services



Life Science Production (LSP) is one of the UK's leading suppliers of off-the-shelf reagents and custom biological products, also offering contract manufacturing services for life science research requirements worldwide.

We provide a first-class service, delivered by knowledgeable, approachable and flexible technical and manufacturing teams, offering full traceability on all products and strong company ethics, ensuring that all clients can be confident that their project is in safe hands.

Some of our quality cell culture products include:

- Foetal Bovine Serum
- Human AB Serum
- Cell Culture Media and Reagents
- CellShip Cell Storage and Transportation Medium
- T-Store Tissue Storage and Transportation Medium


All products are manufactured to a defined specification.





Custom Contract Services

We are a well-known manufacturer and distributor of cell culture reagents. As a leading manufacturer of sera for the past 30 years, this makes us an ideal partner for contract manufacturing and OEM services to the in vitro diagnostic industry.





All our wide range of contract services can be tailored to individual specifications and requirements. In addition, we offer our in-house expertise, two Class 5 ISO 13408-7:2012 cleanrooms and maintain large frozen (-20 °C) and ultra-low temperature (-80 °C) storage facilities. All operations are carried out in our UK facilities, manufacturing with fully disposable, single-use systems.

Examples of custom services include the following:

- Sterile filtration and aseptic filling in a controlled clean room environment
- Flexible packaging in bottles or industrial single-use bags
- Upstream and downstream formulation
- Cell culture media and liquid formulation optimisation
- Research and pilot-scale freezedrying
- Batch sizes from 1 to 1,000 litres

Custom Media Formulation






Custom formulation media manufactured to optimise cell culture manufacturing processes.


Life Science Production (LSP) can manufacture bespoke formulations to exact specifications.

Manufacturing volumes from 10 litres to 1000 litres provide flexibility for media evaluation as well as large-scale production.

LSP has over 30 years of experience in the manufacture of cell culture media and associated reagents.

Biospecimens





Life Science Group: Diverse Biological Specimens and Animal Materials

Life Science Group proudly presents a broad array of biological specimens, constituting one of Europe's largest collections. Our inventory comprises normal and disease-state human donor materials, clinical samples, biofluids, matched sets, and prospective collections. All samples undergo careful collection with the necessary IRB approval and donor consent. With collections in the UK, USA, India, and Europe, our global reach caters to diverse research needs.

We source our animal material from the UK unless specified otherwise. We harvest from animals in fully registered facilities under veterinary supervision, encompassing a wide range of species and strains. Additionally, material is accessible from various agricultural species, expanding the scope of our offerings.

ISO Certified Cleanrooms





Life Science Group has two cleanrooms that guarantee the maintenance of strict sterilised levels, creating controlled manufacturing environments. These spaces provide you with a clean working area tailored for your processes and applications.

Our ISO-classified cleanrooms:

- Adhere to international cleanroom classification standards (ISO 14644)
- Ensure stringent air quality and particulate control
- Comply with cleanroom testing and monitoring standards
- Guarantee strict cleanliness levels for your processes

Bespoke Labelling





Life Science Group possesses the expertise and flexibility to tailor labels according to the specific requirements of our clients.

These labels are meticulously crafted to adhere to all relevant guidelines and specifications, and they can be meticulously designed to accommodate various dimensions and purposes.

The labelled can be edited to have:

- Variable data
- Branded materials
- RFID Technology

Product Sourcing





If an off-the-shelf solution that fits your project requirements is unavailable, we can source and process in bulk material for you. We can pull from a wide range of geographies to suit your needs.

This service is available for both OEM and research-based projects.

Examples of products we can source:

- Blood products
- Serum
- Chick Embryo Extract (CEE)

Antibody Production Services (APS) is one of the UK's leading providers to pharma, bioscience and contract research organisations for flexible and professional contract development of quality antibodies, both monoclonal and polyclonal.

We also provide a wide range of associated services including bulk production, purification, conjugation and optimisation of antibodies, all backed by Life Science Group's commitment to ensuring all ethical requirements are met.

Some of our key services include:

- Custom polyclonal and monoclonal antibodies
- Bulk production
- Purification, conjugation, and optimisation
- Ethically sourced material
- Full technical support and project design
- One of the widest choices of animals for polyclonal antibody production in the UK
- Antibody design for lateral flow
- Antibody sequencing for IP protection

All procedures are performed in secure UK facilities by qualified licensees under veterinary supervision meeting all current UK stringent ethical standards.





Diagnostic Kit Services (DKS) offers tailored, traceable solutions from concept to completion, delivering a variety of high-quality diagnostic kits for both human and veterinary applications. Each kit is fully branded and supported by Life Science Group's commitment to meet every specification comprehensively.

Some of our key services include:

- Complete Flexibility: Our one-stop solution delivers customised kits in any quantity, precisely when needed
- End-to-End Service: From concept to delivery, we ensure consistent quality, traceability, and integrity
- Advanced Manufacturing: Featuring lateral flow technology and scalable production from pilot to full scale
- Specialised Operations: Clean room facilities ensure the highest standards of quality
- Customisation Options: We provide bespoke packaging for medical devices





Contact us:

CliniSciences

info@clinisciences.com www.clinisciences.com



QUALITY | TRACEABILITY | INTEGRITY